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Please read carefully.

EXPLANATION OF SYMBOLS

The following signs and symbols on items, packaging and in the instruction manual bear important information:

This instruction manual belongs to this device. It contains important information about starting up and operation. Read the instruction manual thoroughly. Non-observance of these instructions can result in serious injury or damage to the device.

WARNING These warning notes must be observed to prevent any injury to the user.

CAUTION These notes must be observed to prevent any damage to the device.

NOTE These notes give you useful additional information on the installation or operation.

These test strips for blood glucose tests correspond to the requirements of the EU guidelines 98/79 for in vitro diagnostic devices and bear the CE marking (conformance label) **CE 0483**.

In vitro diagnostic medical device (for external use only)

Discard 6 months after opening Contains sufficient for <n> tests

Do not re-use **LOT** LOT number

Use by **REF** Catalogue number

Storage temperature limitation Manufacturer

Keep away from sunlight

IMPORTANT NOTE Read this leaflet and the instruction manual for the MediTouch® 2 blood glucose monitor carefully before using the MediTouch® test strips. If you still have any questions or require assistance, contact the **MEDISANA** service centre.

INTENDED PURPOSE

The MediTouch® test strips are used with the MediTouch® 2 blood glucose monitor to determine the blood glucose level in fresh capillary blood from the finger tip, alternatively from the ball of the hand or lower arm of adults. Thereby it is a matter of a fast, electro-chemical determination of the blood glucose level. The FAD-binding glucose-dehydrogenase converts the glucose in human blood to gluconolactone. The device measures the current, which is released by this reaction and which is in proportion to the blood glucose volume. The system is intended for external use (in vitro) and can be used for self-testing by persons with diabetes or in clinical settings by healthcare professionals as an aid to monitor the effectiveness of diabetes control.

SAFETY INFORMATION ON STORAGE AND USE

- Do not store the test strips in the refrigerator or freezer.
- Store the test strips in a cool, dry place between 2 °C and 30 °C (35,6 °F - 86 °F).
- Only store the test strips in their original container. Never put them in a new or different container.
- Keep the test strips away from direct sunlight and any other heat source.
- Protect the test strips from high humidity.
- Write down the date you first opened the container on the test strip container.

- The test strips can no longer be used 6 months after opening. Dispose of the remaining test strips in the container together.
- Use each test strip directly after removing it from the container. Close the container immediately after removing a new test strip. Keep the lid sealed at all times.
- A test strip must be used within three minutes of removing it from the container.
- Only touch the test strip with dry, clean hands.
- Protect the test strip from getting contaminated.
- Never bend, cut or modify the test strips.
- Only use MediTouch® 2 test strips.

SAFETY INFORMATION FOR PATIENTS

- Only for use with in vitro diagnostics (for external use only). Do not put in the mouth or swallow.
- Never use the test strips after the expiry date.
- The test strips are only intended to be used once. Never use them twice.
- Medical staff and any others who use this system on more than one patient should be aware that all products or items that come into contact with human blood should be treated as if they are capable of spreading viral diseases, even after cleaning.
- Consult your doctor if you notice symptoms which do not correspond to your test result, despite having followed all the instructions in the manual for the monitor.
- You may only adapt the procedure for using products at home and self-monitoring, if you have first received the appropriate training to do so.

RESTRICTIONS

- MediTouch® blood glucose test strips are for use with a fresh sample of capillary blood. Do not use serum or plasma.
- Do not use the test strips for tests on newborn babies.**
- Extreme humidity can affect the test results. A relative humidity of more than 90 % can lead to incorrect readings.
- The operating temperature of the system is between 10 °C and 40 °C (50 °F - 104 °F). If the temperature is not within this range, incorrect test results can occur.
- Haematocrit: the test results will not be affected by a haematocrit level of between 30% and 55%. A haematocrit level of below 30% can cause an incorrect (too high) test result, a haematocrit level of above 55% can cause an incorrect (too low) test result. If you do not know what your haematocrit level is, consult your doctor.
- Metabolites: the affect of lowering agents depends on the concentration. The following substances may affect the test result - up to the test concentration they will not affect the readings. ➔

| Concentrations of the interferent tested | Glucose Level | 80 mg/dL (4,4 mmol/L) | 250 mg/dL (13,9 mmol/L) | 500 mg/dL (27,8 mmol/L) |
|---|---------------|-----------------------|-------------------------|-------------------------|
| Ascorbic Acid,4 mg/dL(0,26mmol/L) | | 10.89% | -1.76% | 4.55% |
| Ibuprofen,50 mg/dL(2,43mmol/L) | | 3.10% | 2.88% | 4.62% |
| L-Dopa,1,8 mg/dL(0,09mmol/L) | | 10.59% | 7.91% | 4.90% |
| Sodium Salicylate,50 mg/dL(3,12mmol/L) | | -2.59% | 9.42% | -0.84% |
| Tetracycline,1,5 mg/dL(0,03mmol/L) | | -5.32% | 3.81% | 3.20% |
| Tolbutamide,100 mg/dL(3,7mmol/L) | | -2.50% | 12.30% | 0.89% |
| Bilirubin- unconjugated,2,4 mg/dL(0,04mmol/L) | | -2.52% | 4.05% | -0.23% |
| Uric acid,8 mg/dL(0,48mmol/L) | | 2.71% | 9.55% | -1.75% |
| Xylose,4 mg/dL(0,27mmol/L) | | -5.12% | -1.64% | -4.44% |

- Patients undergoing oxygen therapy treatment may receive imprecise test results.

- The test strips can be used at altitudes of up to 3,048 m without having any impact on test results.
- Incorrect test results can occur if the patient is severely dehydrated, suffers from very low blood pressure or is in a state of shock, as well as in cases of hypoglycaemia or hyperosmolarity (with or without ketosis).
- Lipaemia effects: very high levels of blood fat in the blood sample could have a negative impact on certain methodologies. To rule this out, the patient should have his basic level checked in a clinical laboratory test by his doctor before performing the home test. Thereafter, the basic blood glucose level of the patient must be regularly checked and redetermined if necessary.
- Increased cholesterol levels and triglyceride can lead to incorrect test results in certain circumstances.
- Studies have shown that electromagnetism affects the electrical performance and precision of medical equipment and could therefore lead to an incorrect diagnosis.

CONTROL TEST TO ENSURE PRECISE RESULTS

A control test should be performed to ensure precise results of the MediTouch® 2 blood glucose monitor in the following cases:

- You are using the monitor for the first time.
- You have a new container of test strips.
- You have already repeated a test, and the results are still lower or higher than you expected.
- You suspect the monitor or test strips might not be working properly. The monitor has been dropped.
- The test strips have been exposed to a temperature which does not lie within the specified range for the storage conditions (2 °C - 30 °C or 35,6 °F - 86 °F).
- Perform a control test at least once a week.

A control test with control solution is performed to check the precision of the monitor and test strips. Proceed as described in the instruction manual of the MediTouch® 2 blood glucose monitor. The MediTouch® control solution is only to be used with the MediTouch® 2 blood glucose monitor from **MEDISANA**. Other control solutions may cause incorrect test results. If the result of a test with control solution lies within the range printed on the test strip container,

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